

Material

Printer\*



Printer*	Base	Denture 3D+	C&B MFH	Try-In	Model	Model Ortho	Ortho Flex	Ortho Rigid	Ortho IBT	SG (surgical guide)	Cast	Gingiva Mask	Cast	Model 2.0	Tray	Tray	NextDent
NextDent 5100	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
NextDent LCD1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta SOL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta SOL Plus			✓											✓			
Asiga UV Max	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga 4K Pro	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Asiga Ultra	✓		✓	✓	✓	✓		✓	✓	✓	✓				✓	✓	✓
Phrozen CS+			✓							✓							
Phrozen XL 4K		✓	✓					✓	✓	✓				✓			✓
Phrozen Lumii		✓	✓		✓	✓		✓		✓							
Formlabs Form4B					✓					✓					✓	✓	✓
Prusa Medical One		✓	✓				✓			✓				✓	✓	✓	✓
Shining3D Accufab L4D	✓	✓	✓		✓		✓										✓
Shining3D Accufab CEL	✓		✓		✓	✓	✓										✓
Shining3D Accufab D1s		✓	✓		✓	✓	✓										✓
Aextra3D Lumia X1	✓				✓												
DMG DM Desk MC-5			✓	✓													✓
RayShape Edge E2	✓			✓	✓	✓	✓		✓	✓	✓	✓			✓	✓	✓